ONE POT EXTRACTION METHOD FOR MULTIRESIDUE AND POLAR PESTICIDES ANALYSIS: GREEN CHEMISTRY APPROACH CUTTING TIME AND COSTS

S.H. Monteiro¹, C.M. Barbosa¹, B.S.D. Delmondes¹, J.C.C. Silva², J.A. Gonçalves³, G.C.R.M. Andrade¹

¹Biological Institute, Pesticide Residue Laboratory, Reference Laboratory for Research and Analysis of Contaminants in Food and Environment, São Paulo, Brazil, ²Thermo Fisher Scientific Brasil – Nova analítica Imp Exp Ltda, Diadema, Brazil, ³Ministry of Agriculture and Livestock, Residues and Contaminants Laboratory – RCA, Federal Laboratory of Animal and Plant Health and Inspection - LFDA/MG, Pedro Leopoldo, Brazil

Preferred Presentation Type: Oral

Abstract text

Text: Sample preparation methods can pose challenges depending on the physicochemical properties of the pesticides under analysis. Polar compounds, exhibiting higher polarity compared to those detected in multiresidue methods, present unique hurdles such as lower sensitivity, matrix interferents, and the requirement of chromatographic columns with diverse retention mechanisms for detection when it does not use the derivatization approach. In this study, we successfully extracted polar pesticides alongside those typically analyzed in multiresidue methods using a single extraction solvent of acetonitrile and water. For polar compounds like glyphosate, glufosinate, AMPA, fosetyl-Al, maleic hydrazine, and ethephon, a 500 μL aliquot was filtered and analyzed using an Acclaim Trinity Q1 column (2.1 x 100 mm, 3 μm). Subsequently, QuEChERS EN 15662 extraction salts were added to the same centrifuge tube for partitioning and cleaned up with PSA. This facilitated the preparation of samples for determination of other compounds by GC-MS/MS and UPLC-MS/MS. Liquid chromatography coupled with tandem mass spectrometry was employed for detecting both polar and multiresidue pesticides. Although extra caution was needed, such as periodic cleaning runs for the column. The polar method was validated with accuracy ranging between 86% and 129%, and coefficients of variation (CV%) below 19% interday at 0.05 mg/kg, and between 67% and 111% (CV% 18) interday at 0.2 mg/kg. Verification of an additional 273 pesticides yielded accuracy rates of 52% to 137% (CV% 20.3) at 0.02 mg/kg. Sharing the extraction method led to reduced preparation time, solvent consumption, and streamlined laboratory workflows. Overall, the results underscore the efficacy of this extraction method, offering a faster and more cost-effective alternative for integrating polar compounds into laboratory routines.

General

- 1. I confirm that I previewed this abstract and that all information is correct. I accept that the content of this abstract cannot be modified or corrected after final submission and I am aware that it will be published as submitted: Yes
- 2. Submission of the abstract constitutes the consent of all authors to publication (e.g. Workshop website, programmes, other promotions, etc.): Yes
- 3. The Author (person submitting the abstract and presenting at the workshop) warrants and represents that no part of the information and content provided by him/her (Hereafter: the "Content") to EPRW 2024 and INTERPLAN Congress, Meeting & Event Management AG (Hereafter: the "Organisers"), nor the publication of any such Content by the Organisers, on the internet or otherwise, infringes any third-party rights, including but not limited to privacy rights and/or intellectual property rights.: Yes
- 4. The Author grants the Organisers a copyright license to reproduce, publish, translate, distribute, and display the text of the Content on a royalty-free, perpetual, irrevocable nonexclusive basis: Yes
- 5. I herewith confirm that the contact details saved in this system are those of the Author, who will be notified about the status of the abstract. The Author is responsible for informing the other authors (co-authors) about the status of the abstract: Yes
- 6. I understand that the Author making a presentation (oral or poster) must be a registered participant Yes
- 7. I understand that I must select a specific Topic for my abstract allocation. Although the Committee will work hard to honour this selection, this cannot be guaranteed. The Committee reserves the right to change the Topic under which the abstract was originally submitted: Yes
- 8. I herewith confirm that the abstract has not been previously published or presented at another international scientific meeting: Yes